

Product Datasheet

TfR1 (5C9) antibody (orb651364)

Description	Mouse monoclonal antibody to TfR1 (5C9)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide corresponding to residues in Human CD71.
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	~ 120 kDa
Uniprot ID	P02786
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.