

Product Datasheet

GFP antibody (orb651431)

Description	Mouse monoclonal antibody to GFP
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IP, WB
Immunogen	recombinant protein
Preservatives	Mouse IgG1, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	26 kDa
Uniprot ID	P42212
Dilution Range	WB:1:2, 000ICC:1:200
Expiration Date	12 months from date of receipt.