

Product Datasheet

MAP2 antibody (orb652061)

Description	Mouse monoclonal antibody to MAP2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Recombinant Protein
Target	MAP2
Form/Appearance	Mouse monoclonal antibody supplied in crude ascites with 0.09% (W/V) sodium azide.
Note	For research use only
Isotype	IgG1,K
Clonality	Monoclonal
MW	199526
Uniprot ID	P11137
NCBI	NP_114035.2 , NP_001034627.1 , NP_114033.2 , NP_002365.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.