

Product Datasheet

TNFRSF8 antibody (orb652149)

Description	Mouse monoclonal antibody to TNFRSF8
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TNFRSF8
Form/Appearance	Mouse monoclonal antibody supplied in crude ascites with 0.09% (W/V) sodium azide.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	63747
Uniprot ID	P28908
NCBI	NP_001268359.1, NP_001234.2
Dilution Range	WB: 1:100-1000
Expiration Date	12 months from date of receipt.