



Product Datasheet

BHLHE41 antibody (orb652207)


Description	Rabbit polyclonal antibody to BHLHE41
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Synthetic Peptide
Target	BHLHE41
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	50498
Uniprot ID	Q9C0J9
NCBI	NP_110389.1
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.



BHLHE41 antibody (N-term) (Cat. #orb6522...



Fluorescent confocal image of MCF-7 cell...



BHLHE41 antibody (N-term) (Cat. #orb6522...