

Product Datasheet

MKI67 antibody (orb652422)

Description	Mouse monoclonal antibody to MKI67
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthesized Peptide
Target	MKI67
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	358694
Uniprot ID	P46013
Dilution Range	IHC-P: 1:500
Expiration Date	12 months from date of receipt.