

Product Datasheet

KCNT1 antibody (orb652698)

Description Rabbit polyclonal antibody to KCNT1

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC, WB


Note For research use only

Clonality Polyclonal


Uniprot ID **Q5JUK3**

Dilution Range IHC-P: 1:50-100

Expiration Date 12 months from date of receipt.



Western blot analysis of lysate from hum...



Immunohistochemical analysis of paraffin...