

Product Datasheet

p44/42 MAPK (Erk1/2) (Phospho-Thr202/Tyr204) antibody (orb653581)

Description	Mouse monoclonal antibody to p44/42 MAPK (Erk1/2) (Phospho-Thr202/Tyr204)
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Note	For research use only
Clonality	Monoclonal
Dilution Range	WB: 1:500
Expiration Date	12 months from date of receipt.