

Product Datasheet

S100a10 antibody (orb653589)

Description Mouse monoclonal antibody to S100a10

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IF, IP, WB


Note For research use only

Clonality Monoclonal


Uniprot ID **P08207**

Dilution Range ICC: 1:200, IP: 1:100


Expiration Date 12 months from date of receipt.



Western blot detection of S100A10 in C6,...



Immunocytochemistry staining of HeLa cel...



Immunoprecipitation analysis of HeLa cel...