

## Product Datasheet

### LC3A/B antibody (orb653597)

|                            |                                     |
|----------------------------|-------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to LC3A/B |
| <b>Species/Host</b>        | Mouse                               |
| <b>Reactivity</b>          | Human, Rat                          |
| <b>Conjugation</b>         | Unconjugated                        |
| <b>Tested Applications</b> | WB                                  |
| <b>Note</b>                | For research use only               |
| <b>Clonality</b>           | Monoclonal                          |
| <b>Uniprot ID</b>          | <b>Q9H492, Q9GZQ8</b>               |
| <b>Dilution Range</b>      | WB: 1:1000                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.     |