

Product Datasheet

Mouse NOTCH1 Protein (orb653659)

Description	NOTCH1 Interacts with DNER, DTX1, DTX2 and RBPJ/RBPSUH. Also interacts with MAML1, MAML2 and MAML3 w...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	NOTCH1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 55.4 kDa. The protein migrates as 58-68 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_032740.3
Purity	95%
Source	Mouse NOTCH1, His Tag (orb653659) is expressed from human 293 cells (HEK293). It contains AA Ala 19 - Gln 526 (Accession # Q01705-1).
MW	55.4 kDa
NCBI	NP_032740.3
Expiration Date	6 months from date of receipt.