


Product Datasheet

COVID-19 S1 protein antibody (orb653663)

Description	Recombinant human monoclonal antibody recognizes COVID-19 S1 protein antibody.
Species/Host	Human
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Target	S
Preservatives	1 mg/ml. IgG1. Liquid in 0.01M Phosphate Buffered Saline, pH 7.4.
Form/Appearance	1 mg/ml. IgG1. Liquid in 0.01M Phosphate Buffered Saline, pH 7.4.
Concentration	1 mg/ml
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Human Cells
Uniprot ID	P0DTC2
Entrez	43740568
Dilution Range	ELISA: Use at an assay dependent dilution.
Expiration Date	12 months from date of receipt.



Indirect
ELISA
analysis of
COVID-19
S1 P...