

Product Datasheet

COVID-19 Nucleocapsid Protein scFv Antibody (orb653665)

Description Recombinant monoclonal antibody is produced by E.coli expression system with a 6His, Flag tag at the...

Species/Host E. coli

Reactivity Virus

Conjugation Unconjugated

Tested Applications ELISA

Target N

Preservatives 1 mg/ml. IgG1. Liquid in 0.02M Phosphate Buffered Saline, pH 7.4.

Form/Appearance 1 mg/ml. IgG1. Liquid in 0.02M Phosphate Buffered Saline, pH 7.4.

Concentration 1 mg/ml

Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Clonality Monoclonal


Source E.coli

Uniprot ID [P0DTC9](#)

Entrez [43740575](#)

Dilution Range ELISA: Use at an assay dependent dilution.

Expiration Date 12 months from date of receipt.



Indirect
ELISA
analysis of
COVID-19
Nucl...