

Product Datasheet

BID antibody (orb653670)

Description	Rabbit polyclonal antibody to BID
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Recombinant protein of human BID.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.4.
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 25 kDa
Uniprot ID	P55957
Dilution Range	WB: 1:2000~1:5000
Expiration Date	12 months from date of receipt.