

Product Datasheet

ATP1A2 antibody (orb653690)

Description	Rabbit polyclonal antibody to ATP1A2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide, corresponding to Human ATP1A2 .
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2.
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 95kDa
Uniprot ID	P50993
Dilution Range	WB: 1:500~1:1000
Expiration Date	12 months from date of receipt.