

Product Datasheet

NeuN antibody (orb653694)

Description	Rabbit polyclonal antibody to NeuN
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant full length Human NeuN.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1 mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 50 kDa
Uniprot ID	A6NFN3
Dilution Range	WB: 1:500~1:2000 IHC: 1:50~1:200 IF: 1:50~1:200
Expiration Date	12 months from date of receipt.