



Product Datasheet

NMDAR2A (phospho-Tyr1325) antibody (orb6537)

Description	NMDAR2A (phospho-Tyr1325) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human NMDAR2A around the phosphorylation site of Tyr1325 (NF(p-Y)GS)
Target	Phospho-NMDAR2A (Tyr1325)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	165 kDa
Uniprot ID	Q12879
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
mouse
cerebellum
using...



IHC-P
analysis of
rat
cerebellum
using P...