


Product Datasheet

NMDAR2B (phospho-Tyr1070) antibody (orb6538)

Description	NMDAR2B (phospho-Tyr1070) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human NMDAR2B around the phosphorylation site of Tyr1070 (VT(p-Y)GN)
Target	Phospho-NMDAR2B (Tyr1070)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	164 kDa
CAS Number	Glutamate receptor ionotropic, NMDA 2B
Uniprot ID	Q13224
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...