

Product Datasheet

GST Recombinant Protein (orb653871)

Description	Glutathione S- Transferase
Tested Applications	SDS-PAGE
Preservatives	Phosphate buffer (0.1 M, pH 7.4); no azide no glycerol
Concentration	1 mg/ml
Storage	Store at -80°C for long-term storage.
Note	For research use only
Application notes	The protein solution should be gently mixed before use. Avoid freeze/thaw cycles.
Purity	> 95% by SDS-PAGE
Uniprot ID	P08515
Expiration Date	6 months from date of receipt.