

Product Datasheet

NMDAR2B (phospho-Tyr1336) antibody (orb6541)

Description	NMDAR2B (phospho-Tyr1336) antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human NMDAR2B around the phosphorylation site of Tyr1336 (SP(p-Y)AH)
Target	Phospho-NMDAR2B (Tyr1336)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	164 kDa
CAS Number	Glutamate receptor ionotropic, NMDA 2B
Uniprot ID	Q13224
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.