

Product Datasheet

RENT1/hUPF1 Antibody (monoclonal, 11E7) (orb654271)

Description

RENT1/hUPF1 Antibody (monoclonal, 11E7)

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IF, IHC, WB

Applications
Immunogen

A synthetic peptide corresponding to a sequence in the middle region of human RENT1/hUPF1 (578-614aa NMDSMPQLQKLQKDETGELSSADEKRYRALKRT AE), identical to the related mouse and rat sequences.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Mouse IgG2b

Clonality

Monoclonal

MW

130 kDa

Uniprot ID


[Q92900](#)

Dilution Range


Western blot, 0.1-0.5µg/ml, Human, Mouse, Rat
Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, Mouse, Rat
Immunocytochemistry/Immunofluorescence, 2µg/ml, Human

Expiration Date


12 months from date of receipt.



IHC analysis of RENT1/hUPF1 using anti-R...



IHC analysis of RENT1/hUPF1 using anti-R...



IHC analysis of RENT1/hUPF1 using anti-R...