


Product Datasheet

BubR1/BUB1B Antibody (monoclonal, 8B3) (orb654277)


Description

BubR1/BUB1B Antibody (monoclonal, 8B3)


Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	E.coli-derived human BubR1/BUB1B recombinant protein (Position: K26-E448).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Mouse IgG2b
Clonality	Monoclonal
MW	130 kDa
Uniprot ID	O60566
Dilution Range	Western blot, 0.1-0.5µg/ml, Human Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, Rat Flow Cytometry, 1-3µg/1x10 ⁶ cells, Human
Expiration Date	12 months from date of receipt.



Western blot analysis of BubR1/BUB1B using anti-BubR1/BUB1B antibody.



IHC analysis of BubR1/BUB1B using anti-BubR1/BUB1B antibody.



IHC analysis of BubR1/BUB1B using anti-BubR1/BUB1B antibody.