


Product Datasheet


ALKBH5 Antibody (orb654309)

Description	ALKBH5 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
Immunogen	E.coli-derived human ALKBH5 recombinant protein (Position: E81-R180).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	52 kDa
Uniprot ID	Q6P6C2
Dilution Range	Western blot, 0.1-0.25µg/ml, Human Immunocytochemistry/Immunofluorescence, 5µg/ml, Human Direct ELISA, 0.1-0.5µg/ml, Human
Expiration Date	12 months from date of receipt.



IF analysis of ALKBH5 using anti-ALKBH5

...



Western blot analysis of ALKBH5 using an...