

Product Datasheet

CD239/BCAM Antibody (orb654311)

Description

CD239/BCAM Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FC, ICC, IF, IHC, WB

Applications
Immunogen

A synthetic peptide corresponding to a sequence of human CD239/BCAM (LLSLTSTLYLRLRKDDRDASFHCAAHYSLPEGRH).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

70 kDa

Uniprot ID
P50895
Sensitivity


15pg/ml

Dilution Range


Western blot, 0.25-0.5µg/ml, Human, Mouse, Rat Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human Immunocytochemistry/Immunofluorescence, 2µg/ml, Human Flow Cytometry, 1-3µg/1x10⁶ cells, Human

Expiration Date


12 months from date of receipt.



Western blot analysis of CD239/BCAM using anti-CD239/BCAM antibody.



IHC analysis of CD239/BCAM using anti-CD239/BCAM antibody.



IF analysis of CD239/BCAM using anti-CD239/BCAM antibody.