

Product Datasheet

Cd79a Antibody (orb654317)

Description
Cd79a Antibody
Species/Host

Rabbit

Reactivity

Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, FC, IF, IHC, WB

Applications
Immunogen

E.coli-derived mouse Cd79a recombinant protein (Position: L29-P220).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW


54 kDa

Uniprot ID
P11911
Dilution Range


Western blot, 0.25-0.5µg/ml, Mouse, Rat
Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Mouse, Rat
Immunofluorescence, 2µg/ml, Mouse, Rat
Flow Cytometry, 1-3µg/1x10⁶ cells, Mouse
Direct ELISA, 0.1-0.5µg/ml, Mouse

Expiration Date


12 months from date of receipt.



Western blot analysis of Cd79a using anti-Cd79a antibody.



IHC analysis of Cd79a using anti-Cd79a antibody.



IHC analysis of Cd79a using anti-Cd79a antibody.