

Product Datasheet

LDL Receptor/LDLR Antibody (orb654336)

Description

LDL Receptor/LDLR Antibody

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

ELISA, FC, IHC, WB

Applications

Immunogen

E.coli-derived human LDL Receptor/LDLR recombinant protein (Position: Q35-D843).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

130 kDa

Uniprot ID


P01130

Dilution Range


Western blot, 0.25-0.5µg/ml, Human, Rat Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human Flow Cytometry, 1-3µg/1x10⁶ cells, Human, Rat Direct ELISA, 0.1-0.5µg/ml, Human

Expiration Date


12 months from date of receipt.



Western blot analysis of LDLR using anti...



IHC analysis of LDLR using anti-LDLR ant...



IHC analysis of LDLR using anti-LDLR ant...