

Product Datasheet

MBNL1 Antibody (orb654342)

Description
MBNL1 Antibody
Species/Host

Rabbit

Reactivity

Human, Monkey, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, ICC, IF, IHC, WB

Applications
Immunogen

E.coli-derived human MBNL1 recombinant protein (Position: Q89-M388).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW


40 kDa

Uniprot ID
Q9NR56
Dilution Range


Western blot, 0.25-0.5µg/ml, Human, Monkey, Rat
Immunohistochemistry (Paraffin-embedded Section), 2-5µg/ml, Human, Mouse, Rat
Immunocytochemistry/Immunofluorescence, 5µg/ml, Human
Immunofluorescence, 5µg/ml, Human Direct ELISA, 0.1-0.5µg/ml, Human

Expiration Date


12 months from date of receipt.



Western blot analysis of MBNL1 using ant...



IHC analysis of MBNL1 using anti-MBNL1 a...



IHC analysis of MBNL1 using anti-MBNL1 a...