

Product Datasheet

E3 SUMO-protein ligase PIAS1 PIAS1 Antibody (orb654360)

Description

E3 SUMO-protein ligase PIAS1 PIAS1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, FC, ICC, IF, WB

Applications
Immunogen

E.coli-derived human PIAS1 recombinant protein
(Position: Q96-K300).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration
of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After
reconstitution, at 4°C for one month. It can also be
aliquotted and stored frozen at -20°C for six months.
Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive
results. Other applications have not been tested. Optimal
dilutions should be determined by end users. Add 0.2ml
of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW


72 kDa

Uniprot ID
075925
Dilution Range


Western blot, 0.25-0.5µg/ml, Human
Immunocytochemistry/Immunofluorescence, 2µg/ml,
Human Flow Cytometry, 1-3µg/1x10⁶ cells, Human Direct
ELISA, 0.1-0.5µg/ml, Human

Expiration Date


12 months from date of receipt.



Western blot
analysis of
PIAS1 using
ant...



IF analysis of
PIAS1 using
anti-PIAS1
an...



Flow
Cytometry
analysis of
PC-3 cells
us...