

Product Datasheet

TRIM25/EFP Antibody (orb654376)

Description

TRIM25/EFP Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

ELISA, FC, ICC, IF, IHC, WB

Immunogen

E.coli-derived human TRIM25/EFP recombinant protein (Position: M1-A610).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

71 kDa

Uniprot ID

[Q14258](#)

Sensitivity

> 5000 cells


Dilution Range

Western blot, 0.1-0.25µg/ml, Human Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human Immunocytochemistry/Immunofluorescence, 2µg/ml, Human Flow Cytometry, 1-3µg/1x10⁶ cells, Human Direct ELISA, 0.1-0.5µg/ml, Human


Expiration Date

12 months from date of receipt.


Western blot analysis of TRIM25 using anti-TRIM25 antibody



Western blot analysis of TRIM25 using anti-TRIM25 antibody



IHC analysis of TRIM25 using anti-TRIM25 antibody



IHC analysis of TRIM25 using anti-TRIM25 antibody