

Product Datasheet

TIF1 gamma/TRIM33 Antibody (monoclonal, 8I8) (orb654390)

Description

TIF1 gamma/TRIM33 Antibody (monoclonal, 818)

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

FC, ICC, IF, IHC, WB

Immunogen

E.coli-derived human TIF1 gamma recombinant protein (Position: M1001-K1127). Human TIF1 gamma shares 96.1% amino acid (aa) sequence identity with mouse TIF1 gamma.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Mouse IgG2b

Clonality

Monoclonal

MW

150 kDa

Uniprot ID


[Q9UPN9](#)

Dilution Range


Western blot, 0.1-0.5µg/ml, HumanImmunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, HumanImmunocytochemistry/Immunofluorescence, 4µg/ml, HumanFlow Cytometry, 1-3µg/1x10⁶ cells, Human

Expiration Date


12 months from date of receipt.



Western blot analysis of TIF1 gamma/TRIM...



IHC analysis of TIF1 gamma/TRIM33 using ...



IHC analysis of TIF1 gamma/TRIM33 using ...