

Product Datasheet

CTHRC1 Antibody (orb654415)

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat

Tested Applications	ELISA, FC, ICC, IF, IHC, WB
----------------------------	-----------------------------

Form/Appearance	Lyophilized
------------------------	-------------

Storage Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note For research use only

Application notes Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype Rabbit IgG


Clonality	Polyclonal
------------------	------------

MW 26 kDa


Uniprot ID **Q96CG8**

Dilution Range	Western blot, 0.25-0.5µg/ml, Human, Mouse, Rat Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human Immunocytochemistry/Immunofluorescence, 4µg/ml, Mouse Flow Cytometry, 1-3µg/1x10 ⁶ cells, Human Direct ELISA, 0.1-0.5µg/ml, Human
-----------------------	--


Expiration Date	12 months from date of receipt.
------------------------	---------------------------------



Western blot
analysis of
CTHRC1 using
an...



IHC analysis of CTHRC1 using anti-CTHRC1...



Flow
Cytometry
analysis of
HepG2 cells
U...