

## Product Datasheet

### Endothelin A Receptor/ET-A/EDNRA Antibody (orb654421)

## Description

Endothelin A Receptor/ET-A/EDNRA Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

ELISA, FC, IHC, WB

### Applications

### Immunogen

E.coli-derived human Endothelin A Receptor/ET-A/EDNRA recombinant protein (Position: N32-Q381).

### Form/Appearance

Lyophilized

### Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

### Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

### Note

For research use only

### Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

### Isotype

Rabbit IgG

### Clonality

Polyclonal

### MW

50 kDa

### Uniprot ID


**P25101**

### Dilution Range


Western blot, 0.25-0.5µg/ml, Human Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human Flow Cytometry, 1-3µg/1x10<sup>6</sup> cells, Human Direct ELISA, 0.1-0.5µg/ml, Human

### Expiration Date


12 months from date of receipt.



Western blot analysis of Endothelin A Re...



IHC analysis of Endothelin A Receptor/ET...



Flow Cytometry analysis of A549 cells us...