

Product Datasheet

HRP Conjugated Donkey Anti-Rat IgG SABC ELISA Kit [Out of Stock] [Out of stock] (orb654487)

Description	HRP Conjugated Donkey Anti-Rat IgG SABC ELISA Kit
Tested Applications	ELISA
Storage	4°C for one year. Avoid freezing.
Note	For research use only
Application notes	Rabbit IgG refers to the animal origin of the primary antibody, not the origin of the specimen. This kit must be used on primary antibodies from rabbit.
Expiration Date	Please enquire.