

Product Datasheet

HRP Conjugated Donkey Anti-Mouse IgM SABC ELISA Kit [Out of Stock] [Out of stock] (orb654488)

Description	HRP Conjugated Donkey Anti-Mouse IgM SABC ELISA Kit
Tested Applications	ELISA
Storage	4°C for one year. Avoid freezing.
Note	For research use only
Application notes	Mouse IgM refers to the animal origin of the primary antibody, not the origin of the specimen. This kit must be used on primary antibodies from mouse.
Expiration Date	Please enquire.