

Product Datasheet

2019-nCoV 3CLpro, Active protein [Out of stock] (orb654809)

Description Recombinant full-length 2019-nCoV 3CLpro, part of a large replicase polyprotein 1ab (S3264-Q3569), w...

Reactivity Virus

Conjugation Unconjugated

Tested Applications WB

Target 2019-nCoV 3CLpro

Preservatives 50mM sodium phosphate, pH 7.5, 300mM NaCl, 200mM imidazole, 1 mM DTT, 10% glycerol

Storage Store product at -70°C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.

Tag His

Note For research use only

Application notes 2019-nCoV 3CLpro

Source Recombinant full-length 2019-nCoV 3CLpro, part of a large replicase polyprotein 1ab (S3264-Q3569), was expressed in E. coli cells with an N-terminal His tag.

MW 35 kDa

Expiration Date 6 months from date of receipt.