

Product Datasheet

2019-nCoV Nucleocapsid (N) protein [Out of stock] (orb654811)

Description	Recombinant full-length 2019-nCoV Nucleocapsid (N) Protein was expressed in E. coli cells with a C-t...
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Target	2019-nCoV Nucleocapsid (N) Protein
Preservatives	10mM PBS, pH 8.0, 50mM NaCl, 0.05% Na3N.
Storage	Store product at -70°C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.
Tag	His
Note	For research use only
Application notes	2019-nCoV Nucleocapsid (N) Protein
Source	Recombinant full-length 2019-nCoV Nucleocapsid (N) Protein was expressed in E. coli cells with a C-terminal His tag.
MW	52 kDa
Expiration Date	6 months from date of receipt.