

Product Datasheet

2019-nCoV PLpro, Active protein [Out of stock] (orb654812)

Description	Recombinant 2019-nCoV PLpro/papain-like protease, part of a large replicase polyprotein 1ab (E1564-Y...
Reactivity	Virus
Conjugation	Unconjugated
Target	2019-nCoV PLpro
Preservatives	50mM sodium phosphate, pH 7.5, 300mM NaCl, 150mM imidazole, 1mM DTT, 10% glycerol.
Storage	Store product at -70°C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.
Tag	His
Note	For research use only
Application notes	2019-nCoV PLpro
Source	Recombinant 2019-nCoV PLpro/papain-like protease, part of a large replicase polyprotein 1ab (E1564-Y1882), was expressed in E. coli cells with a C-terminal His tag.
MW	38 kDa
Expiration Date	6 months from date of receipt.