

Product Datasheet

Human Novel Coronavirus Nucleoprotein (SARS-CoV-2 N) IgG Antibody ELISA Kit (orb654831)

Description	Human Receptor for advanced glycation end products, Novel Coronavirus Nucleoprotein (SARS-CoV-2 N) I...
Reactivity	Human
Range	Request Information
Target	SARS-CoV-2 N Ab (IgG)
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma
Uniprot ID	P0DTC9
Sensitivity	Request Information
Expiration Date	Please enquire.