

Product Datasheet

Human Novel Coronavirus Nucleoprotein/Spike Glycoprotein S1 (SARS-CoV-2 N/S1) IgG Antibody ELISA Kit (orb654833)

Description	Human Receptor for advanced glycation end products, Novel Coronavirus Nucleoprotein/Spike Glycoprote...
Reactivity	Human
Range	Request Information
Target	SARS-CoV-2 N/S1 Ab (IgG)
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma
Sensitivity	Request Information
Expiration Date	Please enquire.