

Product Datasheet

DNA/ BstEII Digest (orb65515)

Description	The BstEII digest of λ DNA yields the following discrete fragments (in base pairs): 8453*, 7242, 63...
Form/Appearance	liquid (supplied in Jena Bioscience 6x Green Gel Loading Buffer (#PCR-254-gr))
Concentration	200 ng/ μ l
Storage	store at -20 °C
Note	For research use only
Application notes	Preparation: Lambda DNA (48.502 bp) was completely digested by BstEII, phenol extracted, ethanol precipitated, dissolved in 10 mM Tris-HCl (pH 8.0) and 1 mM EDTA. Recommended Load: 5 μ l per lane.
Expiration Date	12 months from date of receipt.