

Product Datasheet

DNA/ HindIII Digest (orb65520)

Description	The HindIII digest of λ DNA yields the following 8 discrete fragments (in base pairs): 23130*, 9416...
Form/Appearance	liquid, DNA fragment mix in 10 mM Tris-HCl pH 7.5, 10 mM EDTA, 8% (w/v) Glycerol, Tartrazine and Xylene Cyanol FF; Color: green
Concentration	100 ng/ μ l
Storage	store at -20 °C
Note	For research use only
Application notes	Preparation: Lambda DNA (48.502 bp) was completely digested by HindIII, phenol extracted, ethanol precipitated, dissolved in 10 mM Tris-HCl (pH 8.0) and 1 mM EDTA. Recommended Load: 5 μ l per lane.
Expiration Date	12 months from date of receipt.