

Product Datasheet

RANK antibody (orb6560)

Description

RANK antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunogen

KLH conjugated synthetic peptide derived from human RANK. Please contact us for the exact immunogen sequence. The peptide is available as orb393922.

Target

RANK

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by Protein A affinity chromatography.

MW

130kDa

Uniprot ID

[Q9Y6Q6](#)

NCBI

[001270951](#), [001257880](#), [11](#)

Entrez


[8792](#)

Dilution Range


IHC-P: 1:100-1:500

Expiration Date


12 months from date of receipt.



IHC-P staining of human breast cancer tissue.



Immunohistochemical staining of paraffin-embedded tissue.



Immunohistochemical staining of paraffin-embedded tissue.