



Product Datasheet

p38 MAPK (phospho-Thr180/Tyr182) antibody (orb6578)


Description	p38 MAPK (phospho-Thr180/Tyr182) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human p38 MAPK around the phosphorylation site of Thr180/Tyr182 (EM(p-T)G(p-Y)VA)
Target	phospho-p38 MAPK (Thr180 + Tyr182)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	41 kDa
Uniprot ID	Q16539
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, ICC=1:100, Flow-Cyt=1µg /test
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human pl...



Western blot analysis of extracts from M...



Flow cytometric analysis of HepG2 Cell u...