

Product Datasheet

Microcystin LR (orb65803)

Description	Cyclic heptapeptide isolated from the freshwater cyanobacteria <i>Microcystis aeruginosa</i> . Potent inhibi...
Conjugation	Unconjugated
Tested Applications	HPLC, NMR
Target	Microcystin LR
Preservatives	Lyophilized solid
Storage	Stable at -25°C for 2 years.
Note	For research use only
Formula	$C_{49}H_{74}N_{10}O_{12}$
Purity	> 95% by HPLC
MW	995.2
CAS Number	[101043-37-2]
Expiration Date	12 months from date of receipt.