

Product Datasheet

Human METAP2 protein (orb6582)

Description	Recombinant of human METAP2 protein (N-6His)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris, 500mM NaCl, 10% glycerol, pH 8.0.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 500mM NaCl, 10% Glycerol, pH 8.0 .
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Protein Sequence	Ala2-Tyr478
Source	Baculovirus
Uniprot ID	P50579
Expiration Date	6 months from date of receipt.