


Product Datasheet

MKK4 (phospho-Ser257/Thr261) antibody (orb6585)

Description	MKK4 (phospho-Ser257/Thr261) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human SEK1/MKK4 around the phosphorylation site of Ser257/Thr261 (VD(p-S)IAK(p-T)RD)
Target	Phospho-MEK4 (Ser257 + Thr261)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	44 kDa
Uniprot ID	P45985
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ce...