

Product Datasheet

MEK5 (phospho-Ser311+Thr315) antibody (orb6586)

Description	MEK5 (phospho-Ser311+Thr315) antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human MEK5 around the phosphorylation site of Ser311+Thr315 (VN(p-S)IAK(p-T)YV)
Target	phospho-MEK5 (Ser311+Thr315)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	49 kDa
Uniprot ID	Q13163
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.