


Product Datasheet

P53 (phospho-Ser46) antibody (orb6598)

Description	P53 (phospho-Ser46) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human P53 around the phosphorylation site of Ser46 (ML(p-S)PD)
Target	phospho-P53 (Ser46)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	43 kDa
Uniprot ID	P04637
Dilution Range	WB=1:500-2000, ICC=1:100, Flow-Cyt=0.2µg /test, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...