



Product Datasheet

P53 (phospho-Thr18) antibody (orb6601)


Description	P53 (phospho-Thr18) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human P53 around the phosphorylation site of Thr18 (QE(p-T)FS)
Target	phospho-P53 (Thr18)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	43 kDa
Uniprot ID	P04637
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, Flow-Cyt=1µg/test, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
mouse
bladder
using ph...



IHC-P
analysis of
Rat bladder
using phos...



IHC-P
analysis of
Human
esophageal
cance...