

Product Datasheet

p70 S6 Kinase Beta (phospho-Thr229) antibody (orb6618)

Description p70 S6 Kinase Beta (phospho-Thr229) antibody

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications FC, WB

Immunogen KLH conjugated synthesised phosphopeptide derived from human RPS6KB2 around the phosphorylation site of Thr228 (TH(p-T)FC)

Target Phospho-p70 S6 Kinase Beta (Thr228)

Preservatives 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance Liquid

Concentration 1mg/ml

Storage Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 59 kDa

Uniprot ID [Q9UBS0](#)

Dilution Range WB=1:500-2000, ELISA=1:5000-10000

Expiration Date 12 months from date of receipt.



Western blot analysis of rat brain lysat...